STANDARD OPERATING PROCEDURE #207 RODENT CASTRATION

- 4.8.6. The fat surrounding the vas deferens and spermatic blood vessels may be gently removed using dry sterile gauze to facilitate cauterization.
- 4.8.7. Cauterize the vas deferens and spermatic blood vessels. Place removed testes aside. Check for bleeding.
- 4.8.8. Use sterile cotton-tipped swabs to gently return any remaining tissues into the scrotum.
- 4.8.9. Repeat with the other testes.
- 4.8.10. Apply topical local analgesic to the incision.
- 4.8.11. Hold the edges of the incision together with forceps and use one drop of tissue glue to close the skin. Alternatively, suture the incision.
- 4.9. Disinfect the instruments between each animal by dipping them in a hot glass bead sterilizer for approximately 30 seconds after removing any blood and debris (let cool completely).
- 4.10. Allow animals to recover in a clean cage. Provide supplemental heat (use a heating disc or pad, heating lamp or incubator) for approximately 30 minutes and monitor the animals until they have fully recovered prior to returning them to their housing room.
- 4.11. Administer analgesics post-operatively as per Rodent Analgesia SOP.

SOP REVISION HISTORY

DATE

NEW VERSION



Rodent Procedure Log



Investigator:	Protocol:						
Procedure: Castration	Performed by:						

Instructions: complete this log for rodent procedures requiring anesthesia, analgesia or post-procedure care (ex. surgeries, experimental infection). Keep the log in the housing room while active and in your files for 3 years for future review by the Quality Assistant and/or the FACC.

active and in your files for 3 years for future review by the Qual	ity Assistant and/or the FACC.					
<u>ANALGESIA</u>	<u>ANESTHESIA</u>	OTHER AGENTS ADMINISTERED				
☐ carprofen: mouse: 20mg/kg, rat: 5-10 mg/kg, SC, every 24 hrs	☐ isoflurane 2-2.5%					
□ buprenorphine: mouse: 0.1mg/kg SC or IP every 4-8 hrs; rat: 0.05mg/kg, SC or IP, every 8-12 hrs	□ ketamine/xylazine/acepromazine*: mouse: 100 mg/kg (K)- 10 mg/kg (X)- 3 mg/kg (A) IP rat: 50 mg/kg (K)- 5 mg/kg (X)- 1 mg/kg (A); IP or IM	□				



Rodent Procedure Log



ANALGESIA

□ carprofen: mouse: 20mg/kg, rat: 5-10 mg/kg, SC, every 24 hrs	
□ buprenorphine: mouse: 0.1mg/kg SC or IP every 4-8 hrs; rat: 0.05mg/kg, SC or IP, every 8-12 hrs	
□ OTHER	Initial the appropriate boxes when complete

	Animal ID Data		Analgesia		SC fluids		Wet food		Time			Remove			
	Animal ID	Date	Day 1	Day 2	Day 3	Day 1	Day 2	Day 3	Day 1	Day 2	Day 3	Day 1	Day 2	Day 3	Sutures (Day 7-10)
1															
2															
3															
4															
5															
6															
7															
8															
9															
10															
11															
12															
13															

Revised: 2014-01-06